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1. Apparatus and reagents

1.1 Apparatus

Name	Producer	Model
Micro CT	NEMO	NMC-200

1.2 Scanning parameters

Project name	Parameters
Radial Field of View (FOV)	100mm
Maximum axial scanning range	250mm
Scanning speed	Max. 4s/bed
Reconstructed pixel size	Min. 2 μm
Spatial resolution	< 7.5 μm @10% MTF

1.3 Softwares

Project name	Software Name
Data acquisition software	Cruiser
Image reconstruction software	Recon
Image display and application software	Avatar

1.4 Major reagents

Name	Producer	Code
Fixative solutions	Servicebio	G1101
Medical alcohol	Hubei Kangjie	/
Sterile gauze	Winner	/

2. Samples

Species	Strain	Model
Quantities	Sample type	Scan date

3. Scanning Procedure

3.1 Scanning: The sample was removed from the fixative, the excess liquid was dried with gauze, the sample was placed in the scanning bed of the instrument and then started to be scanned, the scanning parameter was set to: XXXXXXXXXX, and the original image was obtained after scanning.


3.2 Reconstruction: Reconstruct the original image using the 3D reconstruction software Recon. Adjust the reconstruction parameter to XXXXXXXXX to start the image reconstruction.

3.3 Analysis: Use the data analysis software to analyze the ROI of the target area, select the same area for analysis for all samples, adjust the analysis parameters to: XXXXXXXXXXXXX, get the required parameter values, and export the data.


4. Results

Experimental raw data, including the results of the original graph, analytical data, detailed in the results of the annex, need to be analyzed by the customer and statistical processing.

Example 1: Rat osteoporosis model (femur)

Example 2: Mouse osteoarthritis model (knee joint)

Description	Bone mineral density	Tissue volume	Bone volume	Percent bone volume	Bone surface	Bone surface / volume ratio	Bone surface density	Trabecular pattern factor	Structure model index	Trabecular thickness	Trabecular number	Trabecular separation
Abbreviation	BMD	TV	BV	BV/TV	BS	BS/BV	BS-TV	Tb.Pf	SMI	Tb.Th	Tb.N	Tb.Sp
中文释义	骨矿物质密度	组织体积	骨体积	骨体积分数	骨表面积	骨表面积/体积比	骨表面和组织体积比	骨小梁模式因子	结构模式指数	骨小梁厚度	骨小梁数目	骨小梁分离度
Unit	g/cm ³	mm ³	mm ³	%	mm ²	1/mm	1/mm	1/mm	-	mm	1/mm	mm

For scientific research purposes only, not for clinical diagnosis!