**Servicebio**

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
The screenshot displays a histology software interface for viewing a mouse embryo section. The main window shows a large, detailed cross-section of the embryo, stained with hematoxylin and eosin (H&E). The interface includes a toolbar at the top with various tools for navigation and zooming. A scale bar indicates 2000  $\mu\text{m}$ . A rotation slider is set to -65°. Three red boxes highlight the rotation and flip controls: 'a.任意角度旋转' (arbitrary angle rotation), 'b.水平翻转' (horizontal flip), and 'c.垂直翻转' (vertical flip). The software title is '胎鼠 HE - CaseViewer' and the image is labeled '胎鼠 HE'.

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The screenshot displays the CaseViewer interface with a brain scan image. A central 'Annotation Properties' dialog box is open, and a 'Define Color' dialog box is also visible. Red arrows and numbers 1 through 5 point to specific UI elements. A blue box highlights the 'Show Name' section, and a red box highlights the 'Default Color' section. A blue circle highlights a specific annotation on the brain image.

1. Top right toolbar icon (gear)

2. Top right toolbar icon (gear)

3. Annotation Properties dialog box (gear icon)

4. Left side of the 'Show Name' section in the 'Annotation Properties' dialog box

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Define Color dialog box: 255, 0, 0

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The image shows two software windows. The left window, 'CaseViewer', displays a histological slide with a red circle annotation. A red arrow points to the 'Annotations' toolbar with the instruction: ①使用标注工具，圈出感兴趣的地方，保存 (Use the annotation tool to circle the area of interest and save). The 'Annotations' panel shows 'Annotation 1' with a perimeter of 10.3 mm and an area of 4.7 mm². The right window, '3DHISTECH's Slide Converter', shows the conversion settings. A red box highlights the 'To' dropdown menu with the instruction: ②弹出下图所示窗口，选择保存文件格式 (TIFF、MRXS) (Open the window shown in the figure below, select the save file format (TIFF, MRXS)). The 'Output location' section has 'Local folder' selected, with a red box around the folder path and the instruction: ③选择保存位置 (Select save location). The 'Add slide(s) to queue' button is highlighted with a red box and the instruction: ④添加 (Add). The table below shows the conversion progress:

Input slide	Output slide	Status
MRXS D:\扫描文件\小鼠 HE - Annotation 1	MRXS E:\测试\小鼠 HE - Annotation 1_20200509150820769	Finished

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The screenshot displays the CaseViewer interface with a histology image of a mouse embryo. A red box highlights a region of interest with the following measurements: Perimeter: 8.8 mm, Area: 4.7 mm<sup>2</sup>. A 'Snapshot' dialog box is open, showing options for annotations, scale bar, overview, and slide label. A 'Save As' dialog box is also open, showing file format options.

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**a.显示标注** Show annotations  
**b.显示标尺** Show scale bar  
**c.显示概况** Show slide overview  
**d.显示标签** Show slide label

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
Green Blue Red

The screenshot shows the 'Cell\_line - CaseViewer' software interface. The main window displays a fluorescence microscopy image of cells with blue nuclei, green cytoplasm, and red cell membranes. A 'Pseudo-Color' panel is open on the left, showing a color palette and sliders. A red box highlights the 'Save' button in the Pseudo-Color panel, with the annotation '②选中颜色，保存'. Below the Pseudo-Color panel, a 'Channel' panel is visible, showing three channels: Green, Blue, and Red. A red box highlights the 'Green' channel button, with the annotation '①将光标移动到通道上，右键'. The software interface includes a toolbar at the top with various icons for navigation and zooming, and a scale bar in the top left corner indicating 50  $\mu\text{m}$ .

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